

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

848-2531

LrgHnd Pnch, Cer, 250x3, 90°up, Ej



- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Ejector Pin clears bone chips with each squeeze
- Up 90° Angled Foot Plate
- 250 mm Shaft Length
- 3 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

LrgHndI Punch, Cer, 250x3mm, 90°
up-ang, Ejec

http://buxtonbio.com/lrgHnd_pnchcer250x390upej_1284