

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 848-2562

LrgHnd Pnch, Cer, 250x6, 40°up, Ej



- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Ejector Pin clears bone chips with each squeeze
- Up 40° Angled Foot Plate
- 250 mm Shaft Length
- 6 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

LrgHndI Punch, Cer, 250x6mm, 40°  
up-ang, Ejec

[http://buxtonbio.com/lrgHnd\\_pnchcer250x640upej\\_1238](http://buxtonbio.com/lrgHnd_pnchcer250x640upej_1238)