

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

848-5021

Pnch, Cer, 180x2, 90° up, Ejec, TFP



- Ceramic Black Finish
- 0
- Up 90° Angled Foot Plate
- Foot Plate Thickness Thin
- Chip Handling Ejector
- 180 mm Shaft Length

Sm-Hndl Punch, Cer, 180x2mm, 90° up,
Ejec, TFP

http://buxtonbio.com/pnchcer180x290upejectfp_1291