

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 848-5042

Pnch, Cer, 180x4, 40° up, Ejec, TFP



Photo  
Pending

- Ceramic Black Finish
- 0
- Up 40° Angled Foot Plate
- Cervical; Regular Bite
- Chip Handling Ejector
- 180 mm Shaft Length

Sm-Hndl Punch, Cer, 180x4mm, 40° up,  
Ejec, TFP

[http://buxtonbio.com/pnchcer180x440upejectfp\\_1245](http://buxtonbio.com/pnchcer180x440upejectfp_1245)