

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 848-5452

Pnch, Cer, 230x5, 40° up, Ejec, TFP



- Ceramic Black Finish
- 0
- Up 40° Angled Foot Plate
- Cervical; Regular Bite
- Chip Handling Ejector
- 230 mm Shaft Length

Sm-Hndl Punch, Cer, 230x5mm, 40° up,  
Ejec, TFP

[http://buxtonbio.com/pnchcer230x540upejectfp\\_1247](http://buxtonbio.com/pnchcer230x540upejectfp_1247)