

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 848-8141

Bay-Ker,R-hnd,RR,Cer110x4,90°



- Ceramic
- Original
- Offset Right Handed
- Shaft Length 110mm
- Chip Set Ejector
- 90° Angled Foot Plate

Bay-Ker,R-hnd,RR,Cer110x4,90°

[http://buxtonbio.com/baykerrhndrrcer110x490\\_1718](http://buxtonbio.com/baykerrhndrrcer110x490_1718)