

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 848-8261

Bay-Ker,R-hnd,RR,Cer150x6,90°



- Ceramic
- Original
- Offset Right Handed
- Shaft Length 150mm
- Chip Set Ejector
- Foot Plate Angle u90

Bay-Ker,R-hnd,RR,Cer150x6,90°

[http://buxtonbio.com/baykerrhndrrcer150x690\\_1726](http://buxtonbio.com/baykerrhndrrcer150x690_1726)