

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone:973-560-4848
Fax: 973-560-4999

848-8622

Bay-Ker,L-hnd,RR,Cer150x2,40°



- Ceramic
- Original
- Offset Left Handed
- Shaft Length 150mm
- Chip Set Ejector
- 90° Angled Foot Plate

Bay-Ker,L-hnd,RR,Cer150x2,40°

http://buxtonbio.com/baykerlhndrrcer150x240_1700