

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 86-3542

BAP Punch, 250x4mm, 40° up



- TAP Construction
- Up 40° Angled Foot Plate
- Foot Plate Thickness Regular
- Chip Handling None
- Small Handle
- 250 mm Shaft Length

TAP Sm-Handle Punch, 250x4mm,  
40° up-angled

[http://buxtonbio.com/bap\\_punch\\_250x4mm40up\\_1118](http://buxtonbio.com/bap_punch_250x4mm40up_1118)