

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 86-3952

BAP Punch, 280x5mm, 40° up



- TAP Construction
- Up 40° Angled Foot Plate
- Foot Plate Thickness Regular
- Chip Handling None
- Small Handle
- 280 mm Shaft Length

TAP Sm-Handle Punch, 280x5mm,  
40° up-angled

[http://buxtonbio.com/bap\\_punch\\_280x5mm40up\\_1119](http://buxtonbio.com/bap_punch_280x5mm40up_1119)