

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

86-5011

BAP Punch, 180x1mm, 90°up, TFP



- TAP Construction
- Up 90° Angled Foot Plate
- Foot Plate Thickness Thin
- Chip Handling None
- 180 mm Shaft Length

TAP Sm-Hndl Punch, 180x1mm, 90°
up-angle, TFP

http://buxtonbio.com/bap_punch_180x1mm90uptfp_1192