

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

86-5462

BAP Punch, 230x6mm, 40°up, TFP



- TAP Construction
- Up 40° Angled Foot Plate
- Foot Plate Thickness Thin
- Chip Handling None
- 230 mm Shaft Length

TAP Sm-Hndl Punch, 230x6mm, 40°
up-angle, TFP

http://buxtonbio.com/bap_punch_230x6mm40uptfp_1150