

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 86-8264

BAP PnchR-HdlBay, 150x6, 40°TFP



- Satin
- TAP
- Right Handed Offset
- Shaft Length 150mm
- Chip Set None
- Up 40° Angled Thin Foot Plate

TAP Punch, Rt-Handled Bayonet,  
150x6mm, 40°TFP

[http://buxtonbio.com/bap\\_pnchrhdlbay150x640tfp\\_1624](http://buxtonbio.com/bap_pnchrhdlbay150x640tfp_1624)