

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 863-3261

BAP Punch, 200x6mm, 90°up, Ejec



- TAP Construction
- Up 90° Angled Foot Plate
- Foot Plate Thickness Regular
- Chip Handling Ejector
- Small Handle
- 200 mm Shaft Length

TAP Sm-Hndl Punch, 200x6mm, 90°  
up-ang, Ejec

[http://buxtonbio.com/bap\\_punch\\_200x6mm90upejec\\_1177](http://buxtonbio.com/bap_punch_200x6mm90upejec_1177)