

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

863-3522

BAP Punch, 250x2mm, 40°up, Ejec



- TAP Construction
- Up 40° Angled Foot Plate
- Foot Plate Thickness Regular
- Chip Handling Ejector
- small Handle
- 250 mm Shaft Length

TAP Sm-Hndl Punch, 250x2mm, 40°
up-ang, Ejec

http://buxtonbio.com/bap_punch_250x2mm40upejec_1132