

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

863-3952

BAP Punch, 280x5mm, 40°up, Ejec



- TAP Construction
- Up 40° Angled Foot Plate
- Foot Plate Thickness Regular
- Chip Handling Ejector
- small Handle
- 300 mm Shaft Length

TAP Sm-Hndl Punch, 280x5mm, 40°
up-ang, Ejec

http://buxtonbio.com/bap_punch_280x5mm40upejec_1133