

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 863-5421

BAP Pnch, 230x2, 90° up, Ejec, TFP



- TAP Construction
- Up 90° Angled Foot Plate
- Foot Plate Thickness Thin
- Chip Handling Ejector
- 230 mm Shaft Length

TAP Sm-Hndl Punch, 230x2mm, 90°  
up, Ejec, TFP

[http://buxtonbio.com/bap\\_pnch230x290upejectfp\\_1197](http://buxtonbio.com/bap_pnch230x290upejectfp_1197)