

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

863-7661

TAP, Bay-Kers'n, Up, Ej, 90°, 150x6



Take-Apart Small Kerrison Handle  
Punch Rongeur, Up Bayonet, 150mm  
length, 6.0 bite size, 90° foot angle, Up  
biting, standard foot plate

[http://buxtonbio.com/tapbaykersnupej90150x6\\_1568](http://buxtonbio.com/tapbaykersnupej90150x6_1568)