

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

863-7741

TAP, Bay-Kers'n, Up, Ej, 90°, 180x4



Take-Apart Small Kerrison Handle
Punch Rongeur, Up Bayonet, 180mm
length, 4.0 bite size, 90° foot angle, Up
biting, standard foot plate

http://buxtonbio.com/tapbaykersnupej90180x4_1576