

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

863-8241

TAP, Bay-Ker, R-hnd, Ej, 150x4, 90°



Take-Apart Small Kerrison Handle
Punch Rongeur, Right-Handled
Bayonet, 150mm length, 4.0 bite size,
90° foot angle, Up biting, standard foot
plate

http://buxtonbio.com/tapbaykerrhndej150x490_1625