

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

863-8264

TAP, Bay-Ker, RH, Ej, 150x6, 40° tfp



Take-Apart Small Kerrison Handle  
Punch Rongeur, Right-Handled  
Bayonet, 150mm length, 6.0 bite size,  
40° foot angle, Up biting, thin foot plate

[http://buxtonbio.com/tapbaykerrhej150x640tfp\\_1627](http://buxtonbio.com/tapbaykerrhej150x640tfp_1627)