

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

863-8321

TAP, Bay-Ker, R-hnd, Ej, 180x2, 90°



Take-Apart Small Kerrison Handle  
Punch Rongeur, Right-Handled  
Bayonet, 180mm length, 2.0 bite size,  
90° foot angle, Up biting, standard foot  
plate

[http://buxtonbio.com/tapbaykerrhndej180x290\\_1633](http://buxtonbio.com/tapbaykerrhndej180x290_1633)