

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

867-3941

BAP Punch, Cer, 280x4mm, 90° up

Photo
Pending

- Ceramic Black Finish
- Take Apart Construction
- Up 90° Angled Foot Plate
- Foot Plate Thickness Regular
- Chip Handling None
- Small Handle
- 280 mm Shaft Length

TAP Sm-Hndl Punch, Cer, 280x4mm,
90° up-angle

http://buxtonbio.com/bap_punchcer280x4mm90up_1397