

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 867-5232

BAP Pnch, Cer, 200x3mm, 40° up, TFP



- Ceramic Black Finish
- Take Apart Construction
- Up 40° Angled Foot Plate
- Foot Plate Thickness Thin
- Chip Handling None
- 200 mm Shaft Length

TAP Sm-Hndl Punch, Cer,  
200x3mm, 40° up, TFP

[http://buxtonbio.com/bap\\_pnchcer200x3mm40uptfp\\_1381](http://buxtonbio.com/bap_pnchcer200x3mm40uptfp_1381)