

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 867-7511

BAP Pnch, Cer, U-Bay, 110x1, 90°



Photo  
Pending

- Ceramic
- TAP
- Up Offset
- Shaft Length 110mm
- Chip Set None
- 90° Angled Foot Plate

TAP Sm-Hndl Punch, Cer, Up-Baynt,  
110x1mm, 90°

[http://buxtonbio.com/bap\\_pnchcerubay\\_110x190\\_1747](http://buxtonbio.com/bap_pnchcerubay_110x190_1747)