

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

867-8141

HD, Cer Bay-Ker, R-hdl, 110x4, 90°



- Ceramic
- TAP
- Right Handed Offset
- Shaft Length 110mm
- Chip Set None
- 90° Angled Foot Plate

Bayonet Kerrison rongeur, HD Take-A-Part, Ceramic coated, 110 mm shaft, right-handed, 90° , 4mm bite

http://buxtonbio.com/hdcer_baykerrhdl110x490_1803