

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

867-8244

BAP Pnch, CerR-H-B, 150x4, 40°TFP



- Ceramic
- TAP
- Right Handed Offset
- Shaft Length 150mm
- Chip Set None
- 90° Angled Foot Plate

TAP Punch, CerRt-Handl Baynt,
150x4mm, 40°TFP

http://buxtonbio.com/bap_pnchcerrhb150x440tfp_1813