

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone:973-560-4848  
Fax: 973-560-4999

867-8264

BAP Pnch,CerR-H-B,150x6,40°TFP



- Ceramic
- TAP
- Right Handed Offset
- Shaft Length 150mm
- Chip Set None
- 90° Angled Foot Plate

TAP Punch,CerRt-Handl Baynt,  
150x6mm,40°TFP

[http://buxtonbio.com/bap\\_pnchcerrhb150x640tfp\\_1813](http://buxtonbio.com/bap_pnchcerrhb150x640tfp_1813)