

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 867-8352

BAP Pnch, CerR-Hdl Bay, 180x5, 40°



- Ceramic
- TAP
- Right Handed Offset
- Shaft Length 180mm
- Chip Set None
- 90° Angled Foot Plate

TAP Punch, CerRt-Handled Bayonet,  
180x5mm, 40°

[http://buxtonbio.com/bap\\_pnchcerrhdlbay180x540\\_1821](http://buxtonbio.com/bap_pnchcerrhdlbay180x540_1821)