

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

868-3532

BAP Pnch, Cer, 250x3, 40° up, Ejec



Photo
Pending

- Ceramic Black Finish
- Take Apart Construction
- Up 40° Angled Foot Plate
- Lumbar
- Chip Handling Ejector
- Small Handle
- 250 mm Shaft Length

TAP Sm-Hndl Punch, Cer,
250x3mm, 40° up, Ejec

http://buxtonbio.com/bap_pnchcer250x340upejec_1367