

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 868-3542

BAP Pnch, Cer, 250x4, 40° up, Ejec



- Ceramic Black Finish
- Take Apart Construction
- Up 40° Angled Foot Plate
- Lumbar
- Chip Handling Ejector
- Small Handle
- 250 mm Shaft Length

TAP Sm-Hndl Punch, Cer,  
250x4mm, 40° up, Ejec

[http://buxtonbio.com/bap\\_pnchcer250x440upejec\\_1367](http://buxtonbio.com/bap_pnchcer250x440upejec_1367)