

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

868-4061

BAPF-S Pnch, Cer, 180x6, 90°up, Ej



Photo
Pending

- Ceramic Black Finish
- Take Apart Construction
- Up 90° Angled Foot Plate
- Foot Plate Thickness Regular
- Chip Handling None
- Ferris-Smith Handle
- 180 mm Shaft Length

TAP Frs-Sm Punch, Cer,
180x6mm, 90°up, Ej

http://buxtonbio.com/bapfs_pnchcer180x690upej_1398