

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

868-4451

BAPF-S Pnch, Cer, 230x5, 90°up, Ej



Photo  
Pending

- Ceramic Black Finish
- Take Apart Construction
- Up 90° Angled Foot Plate
- Foot Plate Thickness Regular
- Chip Handling None
- Ferris-Smith Handle
- 230 mm Shaft Length

TAP Frs-Sm Punch, Cer,  
230x5mm, 90°up, Ejec

[http://buxtonbio.com/bapfs\\_pnchcer230x590upej\\_1400](http://buxtonbio.com/bapfs_pnchcer230x590upej_1400)