

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

868-5032

BAPPnch, Cer, 180x3, 40° up, Ej, TFP



- Ceramic Black Finish
- Take Apart Construction
- Up 40° Angled Foot Plate
- Foot Plate Thickness Thin
- Chip Handling Ejector
- 180 mm Shaft Length

TAP Sm-Hdl Punch, Cer, 180x3mm, 40°
up, Ejec, TFP

http://buxtonbio.com/bappnchcer180x340upejtfp_1380