

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

868-5051

BAPPnch, Cer, 180x5, 90° up, Ej, TFP



- Ceramic Black Finish
- Take Apart Construction
- Up 90° Angled Foot Plate
- Foot Plate Thickness Thin
- Ejector
- 180 mm Shaft Length

TAP Sm-Hdl Punch, Cer, 180x5mm, 90°
up, Ejec, TFP

http://buxtonbio.com/bappnchcer180x590upejtfp_1429