

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

868-7561

TAP, Bay-Ker, Up, Ej, Cer90°, 110x6



Take-Apart Small Kerrison Handle ,  
Ceramic Coated Punch Rongeur, Up  
Bayonet, 110mm length, 6.0 bite size,  
90° foot angle, Up biting, standard foot  
plate

[http://buxtonbio.com/tapbaykerupejcer90110x6\\_1750](http://buxtonbio.com/tapbaykerupejcer90110x6_1750)