

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

868-7632

TAP, Bay-Ker, Up, Ej, Cer40°, 150x3



Take-Apart Small Kerrison Handle ,  
Ceramic Coated Punch Rongeur, Up  
Bayonet, 180mm length, 3.0 bite size,  
40° foot angle, Up biting, standard foot  
plate

[http://buxtonbio.com/tapbaykerupejcer40150x3\\_1759](http://buxtonbio.com/tapbaykerupejcer40150x3_1759)