

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

868-77-5162

BAPnch, Cer, U-Bay, 18x2, LCUp, TFP



Photo  
Pending

- Ceramic
- TAP
- 180 mm Shaft Length
- Thin Foot Plate

Take-Apart Small Kerrison Handle ,  
Ceramic Coated Punch Rongeur, Up  
Bayonet, 180mm length, 2.0 bite size,  
40° foot angle, Up biting, thin foot plate

[http://buxtonbio.com/bapnchcerubay18x2lcuptfp\\_1851](http://buxtonbio.com/bapnchcerubay18x2lcuptfp_1851)