

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

868-8152

TAP, Bay-Ker, RH, Ej, Cer110x5, 40°



Photo  
Pending

Take-Apart Small Kerrison Handle ,  
Ceramic Coated Punch Rongeur,  
Right-Handled Bayonet, 110mm  
length, 5.0 bite size, 40° foot angle, Up  
biting, standard foot plate

[http://buxtonbio.com/tapbaykerrhejcer110x540\\_1807](http://buxtonbio.com/tapbaykerrhejcer110x540_1807)