

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

868-8251

TAP, BayKer, R-H, Ej, Cer, 90° 150x5



Take-Apart Small Kerrison Handle ,
Ceramic Coated Punch Rongeur,
Right-Handled Bayonet, 180mm
length, 5.0 bite size, 90° foot angle, Up
biting, standard foot plate

http://buxtonbio.com/tapbaykerrhejcer90150x5_1814