

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

868-8264

TAP, BayKer, RH, EjCr150x6, 40° tfp



Take-Apart Small Kerrison Handle ,
Ceramic Coated Punch Rongeur,
Right-Handled Bayonet, 180mm
length, 6.0 bite size, 40° foot angle, Up
biting, thin foot plate

http://buxtonbio.com/tapbaykerrhejcr150x640tfp_1816