

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

868-8622

TAP, Bay-Ker, LH, Ej, Cer150x2, 40°



Photo  
Pending

Take-Apart Small Kerrison Handle ,  
Ceramic Coated Punch Rongeur, Left-  
Handled Bayonet, 180mm length, 2.0  
bite size, 40° foot angle, Up biting,  
standard foot plate

[http://buxtonbio.com/tapbaykerljejcer150x240\\_1787](http://buxtonbio.com/tapbaykerljejcer150x240_1787)