

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 868-8732

TAP, Bay-Ker, LH, Ej, Cer180x3, 40°



Take-Apart Small Kerrison Handle ,  
Ceramic Coated Punch Rongeur, Left-  
Handled Bayonet, 180mm length, 3.0  
bite size, 40° foot angle, Up biting,  
standard foot plate

[http://buxtonbio.com/tapbaykerljejcer180x340\\_1795](http://buxtonbio.com/tapbaykerljejcer180x340_1795)