



# DAHMEN-AIDED SURGERY

## DAHMEN OSTEOTOMES, 300mm

	<b>STRAIGHT</b>	<b>ANGLED</b>
4 mm	77-1104	77-1504
5 mm	77-1105	77-1505
6 mm	77-1106	77-1506
8 mm	77-1108	77-1508
10 mm	77-1110	77-1510
12 mm	77-1112	77-1512
15 mm	77-1115	77-1515
20 mm	77-1120	77-1520
22 mm	77-1122	77-1522
25 mm	77-1125	77-1525

## DAHMEN GOUGES, 300mm

	<b>STRAIGHT</b>	<b>ANG. BACK</b>	<b>ANG. FWD</b>
4 mm	77-2004	77-2204	77-2504
5 mm	77-2005	77-2205	77-2505
6 mm	77-2006	77-2206	77-2506
8 mm	77-2008	77-2208	77-2508
10 mm	77-2010	77-2210	77-2510
12 mm	77-2012	77-2212	77-2512
15 mm	77-2015	77-2215	77-2515
20 mm	77-2020	77-2220	77-2520
22 mm	77-2022	77-2222	77-2522
25 mm	77-2025	77-2225	77-2525

Enjoy the benefits of slim, lightweight 30 cm long osteotomes with linen phenolic handles as offered by Buxton. Their extended length provides a distinct advantage when treating large patients and for various posterior approach procedures. While originally offered for spinal applications the gouges are sometimes employed to remove cement and debris within the intramedullary canal.

Osteotomes are available straight or with a 45° angled head, while gouges are additionally offered as straight, forward or back angled. To aide in the broadest range of applications, widths are available from 4 mm to 25 mm. Dahmen-Lexer osteotomes and gouges are additionally offered as alternatives to the leptosomic Dahmen. These straight, stouter alternatives come in seven widths ranging from 5 mm to 30 mm.

## DAHMEN-LEXER, 225mm

	<b>OSTEOTOME</b>	<b>GOUGE</b>
5 mm	77-4005	77-4205
7 mm	77-4007	77-4207
10 mm	77-4010	77-4210
15 mm	77-4015	77-4215
20 mm	77-4020	77-4220
25 mm	77-4025	77-4225
30 mm	77-4030	77-4230

For a Catalog featuring our Spine Instruments, call or visit

1-800-427-5351

www.buxtonbio.com

**BUXTON BIOMEDICAL, INC.**

11 Melanie Lane, East Hanover NJ 07936

www.buxtonbio.com • info@buxtonbio.com